

What size is needed?

The cooling compression bandage contracts when filled with water and thus becomes "smaller", especially in width, but also slightly in height.

The right size for your horse: The decisive factor is the distance from the carpal joint to the ground. This is how the size is calculated:


Measure the distance from the bottom edge of the carpal joint to the floor.

The applied bandage extends upwards at the back approximately 3 to 4 cm above the measuring point between the cannon bone and the joint. There should still be at least 1 cm of space above the ground.

This results in:

For bandage size **M medium** the distance from the carpal joint to the ground must be **39 (minimum) to 49 cm**.

For the bandage size **L large** the distance from the carpal joint to the ground must be **at least 50 cm**.



Our two sets: E-cooler box or cooling bucket

The sets contain everything you need for one leg (left or right). Options (e.g. a bandage for the other leg) and spare parts can be ordered separately.

Cooling compression bandage set cooling bucket

The **cooling bucket with hand pump** is particularly suited for quick and frequent use for regeneration after training or competition and for short-term cooling of the horse's leg in case of swelling, strains, bruises and/or inflammation below the carpal joint. Since the water hoses are clicked off after being filled with water, the horse can move freely in the box. **Additional cooling compression bandages for cooling multiple legs can be purchased separately.**



Cooling compression bandage set E-cooler box

The **electric cooler** comes with a powerful power bank. A continuous pump continuously pumps ice water through the cooling compression bandage, creating a soothing massage effect. This method allows for a longer cooling time, although after a certain period of time, some of the water may need to be replaced or new ice added. With this method, the horse must be tied up for the duration of the treatment.



Order form

Name of purchaser

Address

Telephone & email

Signature



Cooling compression bandage

Set electric cool box

Set consists of: E-cooling box, water hoses with suitable connectors, instruction manual, vetal transport bag for e-cooler with wheels, power bank with charging cable, pump with USB cable, Switch with USB ports, water drainage set and a cooling compression bandage of size



			Number of sets
or	medium left	1111M-L	<input type="text"/>
or	medium right	1111M-R	<input type="text"/>
or	large left	1111 L-L	<input type="text"/>
or	large right	1111 L-R	<input type="text"/>

Options electric coolers

	Number
Ceiling belt with guide loops (3 short and 1 long Velcro straps)	0170 <input type="text"/>
Hose insulation optional, 2 x 3 meters	0130 <input type="text"/>

Spare parts for electric coolers

Passive cool box 26 l	0111	<input type="text"/>
Hose with suitable connector, L300, 1 piece	0122	<input type="text"/>
Valve closure with tab, 1 piece	0140	<input type="text"/>
vetal transport bag mobile	0151	<input type="text"/>
Power bank with charging cable	0160	<input type="text"/>
USB water pump	0161	<input type="text"/>
USB 2.0 cable with on/off switch	0162	<input type="text"/>
Water drainage set	0123	<input type="text"/>

Set of cooling bucket

Set consists of: cooling bucket, water supply hose with integrated hand pump, water return hose, instructions for use, vetal transport bag for cooling bucket and one cooling compression bandage of size



			Number of sets
or	medium left	1110M-L	<input type="text"/>
or	medium right	1110M-R	<input type="text"/>
or	large left	1110L-L	<input type="text"/>
or	large right	1110L-R	<input type="text"/>

Spare parts cooling bucket

	Number
Cooling bucket	0110 <input type="text"/>
Water supply hose with integrated hand pump	0120 <input type="text"/>
Water return hose	0121 <input type="text"/>
Valve closure with tab, 1 piece	0140 <input type="text"/>
vetal transport bag for cooling buckets	0150 <input type="text"/>
Hose connection cover with lock nut	0141 <input type="text"/>

Options to match both sets



Cooling compression bandage

		Number
medium left	1120M-L	<input type="text"/>
medium right	1120M-R	<input type="text"/>
large left	1120 L-L	<input type="text"/>
large right	1120 L-R	<input type="text"/>



Insulating bandage

medium	1130M	<input type="text"/>
large	1130L	<input type="text"/>